

MINI & MAXI FLANGE ACCESSORIES **Gaskets**

Nordson MEDICAL's gaskets maintain seal integrity, have complete material traceability, and are proven to have excellent wear performance. Material options include Class VI compliant PTFE/EPDM and platinum-cured silicone in gasket sizes that fit both 1-1/2" (Maxi) and 3/4" (Mini) flanges. Keeping customer requirements in mind, these components were manufactured for applications requiring the highest sanitary standards and chemical purity. Supported with a complete regulatory package, our selection of gaskets deliver on our promise to ensure reliable and consistent quality with our products.

FEATURES & BENEFITS

- Maintains seal integrity
- Class VI-compliant materials
- 100% virgin raw materials & colorants
- Excellent wear performance
- Complete material traceability

SPECIFICATIONS

- **Materials:**
Platinum-cured silicone, Polytetrafluoroethylene (PTFE)/
Ethylene propylene diene monomer (EPDM)
- **Flange Size:**
 - 3/4" (19.0 mm) - Mini
 - 1-1/2" (38.1 mm) - Maxi
 - 3" (76.2 mm)
 - 4" (101.6 mm)
 - 6" (152.4 mm)

GASKETS



GMN75-PX

Molded Gasket for Mini Flange SFMN Series Fittings, Platinum-Cured Silicone



G40MPGR-TF150-E

Molded Tuf-Flex® Gasket for Maxi Flange SFMX Series Fittings, PTFE



GMN450-PX

Molded Gasket for Mini Flange SFMN Series Fittings, Platinum-Cured Silicone




G424TRX-0300

Molded Gasket for 3 inch Flange Fittings, Platinum-Cured Silicone



GMX455-PX

Molded Gasket for Maxi Flange SFMX Series Fittings, Platinum-Cured Silicone




G40RX-PX-0400

Molded Gasket for 4 inch Flange Fittings, Platinum-Cured Silicone



GMX150-PX

Molded Gasket for Maxi Flange SFMX Series Fittings, Platinum-Cured Silicone



GM4749TRX-0600

Molded Gasket for 6 inch Flange Fittings, Platinum-Cured Silicone