





BILL OF MATERIAL:

WHITE POLYESTER: 325M VALOX

BLUE POLYESTER: 325M VALOX WITH CPBT02803
BLUE VALOX COLORANT

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY
OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN
PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY	DMC	DATE	5/14/2013
APPROVED BY		DATE	
MATERIAL			
FINISH			
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		THIRD ANGLE PROJECTION 	

HR® Halkey-Roberts®		
TITLE TUBING CLAMP, LARGE		
SIZE B	DRAWING NO. C340TCLX	REV U
SCALE 1:1	DO NOT SCALE DRAWING	SHEET 2 OF 2