




DRAWING NO. CV243001VG	REV C	Customer Drawing	REVISION HISTORY			
			REV	DESCRIPTION (ECO)	DATE	APPROVED
			B	ADD DIMENSION .280, AND REPOSITION .295 ECO 7579	07/20/17	AK
		C	UPDATED STEM MATERIAL FROM CLS2600HR SILICONE TO LR300360800 SILICONE, ECO 13219	03/18/26	JS	

BODY	STEM	PLUG
2212RHT1118C040 CLEAR PVC GAMMA	LR300360800 BLACK SILICONE K71008 COLORANT	H0500HN CLEAR POLYPROPYLENE

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY	DMC	DATE	1/24/2013	HR® Halkey-Roberts®	
APPROVED BY		DATE			
MATERIAL					
FINISH				TITLE	CATHETER VALVE ASSEMBLY
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		THIRD ANGLE PROJECTION		SIZE	DRAWING NO.
				A	CV243001VG
				SCALE	DO NOT SCALE DRAWING
				1:1	SHEET 2 OF 2