

- 2. LETTERS TO BE 3/32 HIGH & .005 MIN. RAISED, AS SHOWN.
- 3. CRITICAL DIMENSIONS ARE CIRCLED, AND SUBJECT TO INSPECTION.
- 4. DIMENSIONS IN PARENTHESES ARE FOR REFERENCE ONLY.
- 5. ALL DIMENSIONS ARE IN INCHES WITH NOMINAL, DIMENSIONS IN MILLIMETERS BELOW.

O.D. FLASH: 0.060 MAX. (1.5MM)

GATE VESTIGE: 0.100 MAX. (2.5MM)

EMBEDDED PARTICLES, BUBBLES, SPOTS, STAINS: 0.3mm² MAX. NO MORE THAN 5 PER PART.

DE	
■ =•	DRAW

	DRAWN BY	M.S.	10-13-97 DATE	/ [/	Halkey-Roberts	
	APPROVED BY MATERIAL	JW PVC	-	<u> </u>	R Haikey-Roberts	4
MOLDED PARTS: GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA	FINISH			TITLE	FLO - CHECK VALVE	
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE	UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES:		THIRD ANGLE PROJECTION	SIZE	C1020AFX)
PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.	ANGULAR: ± 1° TWO PLACE DECIM THREE PLACE DECIM			SCALE 1:1	DO NOT SCALE DRAWING SHEET 1 OF 1	