

DRAWING NO.

CV15900LDB


REV

A

Customer Drawing

Revision History

REV	DESCRIPTION (ECO)	DATE	APPROVED
A	NEW DRAWING ECN 3843	8/06/08	P.M.
-	-	-	-



DRAWN BY

DMC

DATE

APPROVED BY

P.M.

DATE

MATERIAL

SEE SHEET 2

FINISH


UNLESS OTHERWISE SPECIFIED:

DIMENSIONS ARE IN INCHES

TOLERANCES:

ANGULAR: $\pm 1^\circ$ TWO PLACE DECIMAL $\pm .01$ THREE PLACE DECIMAL $\pm .005$

THIRD ANGLE PROJECTION



Halkey-Roberts

TITLE

HAND PUMP, ASSY.
W/ LUER TAPER END

SIZE

A

SCALE

1:1

DRAWING NO.

CV15900LDB

REV

A

DO NOT SCALE DRAWING

SHEET

1 OF 2

MOLDED PARTS:

GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA



PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWING NO. CV15900LDB	REV A	Customer Drawing	Revision History			
			REV	DESCRIPTION (ECO)	DATE	APPROVED
			A	SEE SHEET ONE	-	-
			-	-	-	-

BILL OF MATERIAL

V243101PCG	V157212DB, V15722DB	V15940LPGC
CATHETER VALVES	HAND PUMP, RIGHT AND LEFT	FLOW CONTROL VALVE W/ CATHETER VALVE
POLYPROPYLENE GAMMA RESISTANT	MEDICAL GRADE PVC, DARK BLUE	POLYPROPYLENE GAMMA RESISTANT

DRAWN BY DMC		DATE		 Halkey-Roberts				
APPROVED BY P.M.		DATE						
MATERIAL AS LISTED		TITLE HAND PUMP, ASSY. W/LUER TAPER END						
FINISH		UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		THIRD ANGLE PROJECTION 		SIZE A	DRAWING NO. CV15900LDB	REV A
MOLDED PARTS: GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA		PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.		SCALE 1		DO NOT SCALE DRAWING		SHEET 2 OF 2