


BILL OF MATERIAL					
SUFFIX "X"	STEM RETAINER	MALE LUER	LUER LOCK BODY	TUBE CONNECTOR	STEM
251001004	P5M6K080 POLYPROPYLENE CLEAR	P5M6K080 POLYPROPYLENE CLEAR	RX1805451118 POLYCARBONATE MAKROLON CLEAR	RX1805 451118 POLYCARBONATE MAKROLON CLEAR	LR3003/40 SILICONE WACKER 40 DUROMETER CLEAR
251001004R	P5M6K080 POLYPROPYLENE CLEAR	P5M6K080 POLYPROPYLENE CLEAR	RX1805451118 POLYCARBONATE MAKROLON W/2% RED COLORANT	RX1805451118 POLYCARBONATE MAKROLON CLEAR	LR3003/40 SILICONE WACKER 40 DUROMETER CLEAR
251001004B	P5M6K080 POLYPROPYLENE CLEAR	P5M6K080 POLYPROPYLENE CLEAR	RX1805451118 POLYCARBONATE MAKROLON W/4% BLUE COLORANT	RX1805451118 POLYCARBONATE MAKROLON CLEAR	LR3003/40 SILICONE WACKER 40 DUROMETER CLEAR
251001004G	P5M6K080 POLYPROPYLENE CLEAR	P5M6K080 POLYPROPYLENE CLEAR	RX1805451118 POLYCARBONATE MAKROLON W/4% GREEN COLORANT	RX1805451118 POLYCARBONATE MAKROLON CLEAR	LR3003/40 SILICONE WACKER 40 DUROMETER CLEAR
251001004Y	P5M6K080 POLYPROPYLENE CLEAR	P5M6K080 POLYPROPYLENE CLEAR	RX1805451118 POLYCARBONATE MAKROLON W/2% YELLOW COLORANT	RX1805 451118 POLYCARBONATE MAKROLON CLEAR	LR3003/40 SILICONE WACKER 40 DUROMETER CLEAR

**MOLDED PARTS:**  
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

**PROPRIETARY AND CONFIDENTIAL**  
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY  
OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY  
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN  
PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY	S.Rogier	DATE	12/2/14
APPROVED BY	DK	DATE	06/12/24
MATERIAL	See Sheet 2.		
FINISH	N/A		
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		THIRD ANGLE PROJECTION	



<h2>HR® Halkey-Roberts®</h2>	
TITLE <b>MALE SHUT OFF VALVE 6.6mm</b>	
SIZE <b>A</b>	DRAWING NO. <b>C251001004X</b> REV <b>D</b>
SCALE <b>1:1</b>	DO NOT SCALE DRAWING
SHEET <b>2 OF 2</b>	