


BILL OF MATERIAL				
STEM RETAINER	MALE LUER	LUER LOCK BODY	LUER CONNECTOR	STEM
CLEAR POLYPROPYLENE P5M6K080	CLEAR POLYPROPYLENE P5M6K080	CLEAR POLYCARBONATE MAKROLON RX 1805 451118	CLEAR POLYCARBONATE MAKROLON RX 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY S.Rogier	DATE 2/5/2015	HR® Halkey-Roberts®
APPROVED BY DK	DATE 6/10/24	
MATERIAL SEE SHEET 2.		
FINISH N/A	UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005	THIRD ANGLE PROJECTION 
SIZE A		TITLE MALE SHUT OFF VALVE FLL
SCALE 1:1		DRAWING NO. C251001002
DO NOT SCALE DRAWING		REV E
SHEET 2 OF 2		