


BILL OF MATERIAL				
STEM RETAINER	MALE LUER	LUER LOCK BODY	BARB CONNECTOR	STEM
CLEAR POLYPROPYLENE P5M6K080	CLEAR POLYPROPYLENE P5M6K080	CLEAR POLYCARBONATE MAKROLON RX 1805 451118	CLEAR POLYCARBONATE MAKROLON RX 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY
OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN
PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY S.Rogier	DATE 2/5/15	HR® Halkey-Roberts®	
APPROVED BY D.K	DATE 06/12/24		
MATERIAL See Sheet 2.			
FINISH N/A	UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005	THIRD ANGLE PROJECTION 	TITLE MALE SHUT OFF VALVE .125 BARB
		SIZE A	DRAWING NO. C251001003
		SCALE 1:1	REV D
		DO NOT SCALE DRAWING	
		SHEET 2 OF 2	