


BILL OF MATERIAL

VALVE BODY: 2212RHT1118C040 PVC CLEAR
BONDABLE ADAPTOR: 2212RHT1118C040 PVC CLEAR
STEM: CLS2600HR SILICONE KIT, W/ K71008 COLORANT BLACK
PLUG: H0500HN POLYPROPYLENE CLEAR

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY	DMC	DATE	7/20/09
APPROVED BY		DATE	
MATERIAL			
FINISH			
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005	THIRD ANGLE PROJECTION		



HR® Halkey-Roberts®	
TITLE LUER LOCK VALVE W/ Bondable Adaptor	
SIZE A	DRAWING NO. C247200001
SCALE 1:1	DO NOT SCALE DRAWING
REV F	
SHEET 2 OF 2	