


SUFFIX	LUER BODY	BARBED BASE	STEM
024	CLEAR POLYCARBONATE MAKROLON RX1805451118	CLEAR POLYCARBONATE MAKROLON RX1805451118	BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DR COLOR K75238

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY
OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN
PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY KJS	DATE 6/4/15	HR® Halkey-Roberts®		
APPROVED BY	DATE			
MATERIAL				
FINISH		TITLE 1/8" BARBED SWABABLE VALVE		
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		SIZE A	DRAWING NO. C245508XXX	REV E
THIRD ANGLE PROJECTION 		SCALE 2:1	DO NOT SCALE DRAWING	
		SHEET 2 OF 2		