

BILL OF MATERIAL:

PART NUMBERS	BODY FEMALE	BODY MALE THREADED	STEM SWABABLE
245205024Q	POLYCARBONATE, YELLOW	POLYCARBONATE CLEAR	SILICONE, BLUE COLORANT
245205024Y	POLYCARBONATE, YELLOW	POLYCARBONATE YELLOW	SILICONE, BLUE COLORANT
245205524Q	POLYCARBONATE, YELLOW	POLYCARBONATE, CLEAR	SILICONE, CLEAR
245205524Y	POLYCARBONATE, YELLOW	POLYCARBONATE YELLOW	SILICONE, CLEAR


DRAWN BYDCDATE03/01/23

APPROVED BYSRDATE03/01/23

MATERIALSEE NOTES

FINISH

UNLESS OTHERWISE SPECIFIED:  
DIMENSIONS ARE IN INCHES  
TOLERANCES:  
ANGULAR: ± 1°  
TWO PLACE DECIMAL ± .01  
THREE PLACE DECIMAL ± .005

THIRD ANGLE PROJECTION

HR® Halkey-Roberts®

TITLESWABABLE VALVE, HD

SIZEB

DRAWING NO.C245205X24X

REVC

SCALE4:1

DO NOT SCALE DRAWING

SHEET2 OF 2

MOLDED PARTS:  
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL  
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY  
OF HALKEY-ROBERTS CORPORATION, ST. PETERSBURG, FL. ANY  
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN  
PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.