

SUFFIX	FEMALE HOUSING	MALE HOUSING	STEM
60	BLUE COPOLYESTER (TRITAN MX731)	CLEAR COPOLYESTER (TRITAN MX731)	SILICONE RUBBER (WALKER LR3003/40, BLUE)

DRAWN BYKS

APPROVED BYAK


MATERIALSEE TABLE

FINISHN/A

DATE3/14/23


DATE3/14/23

UNLESS OTHERWISE SPECIFIED:  
DIMENSIONS ARE IN INCHES [MM]  
TOLERANCES:  
ANGULAR: ± 1°  
TWO PLACE DECIMAL ± .01  
THREE PLACE DECIMAL ± .005

THIRD ANGLE PROJECTION

SIZEB

SCALE4:1

Halkey-Roberts

TITLE  
NEUTRAL VALVE

DRAWING NO. C2452010XX

DO NOT SCALE DRAWING

REV  
B

SHEET  
2 OF 2