



BILL OF MATERIAL:

1ST SUFFIX	2ND SUFFIX	LUER BODY	CONNECTOR	STEM
024	FT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238
024R	FT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118 2% RED CLARIANT TRANSPARENT		
024B	FT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118 4% BLUE CLARIANT TRANSPARENT		
524	FT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40