


Customer Drawing

REV	DESCRIPTION (ECO)	DATE	APPROVED
D		07/20/18	AK
E	UPDATE NOTE 3, & 4, ECO 8301	10/09/18	AK



HR® Halkey-Roberts®			
TITLE			
T-CONNECTOR ASSEMBLY			
SIZE	DRAWING NO.		REV
A	C245510024		E
SCALE	DO NOT SCALE DRAWING	SHEET	1 OF 3
2:1			

SHEET	2 OF 3
-------	--------

DRAWING NO. C245510024	REV E	Customer Drawing	REVISION HISTORY			
			REV	DESCRIPTION (ECO)	DATE	APPROVED
			D		07/20/18	AK
			E	UPDATE NOTE 3, & 4, ECO 8301	10/09/18	AK

BILL OF MATERIAL				
CONFIGURATION	T-PORT BASE	SWABABLE CAP	SWIVEL LOCK	STEM
245510024	CLEAR POLYCARBONATE MAKROLON RX 1805 451118	CLEAR POLYCARBONATE MAKROLON RX 1805 451118	POLYCARBONATE CALIBRE MEGARAD 2081-15	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY	SR	DATE	12/9/15	HR® Halkey-Roberts®	
APPROVED BY		DATE			
MATERIAL					
FINISH				TITLE	T-CONNECTOR ASSEMBLY
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		THIRD ANGLE PROJECTION		SIZE	DRAWING NO.
				A	C245510024
				SCALE	DO NOT SCALE DRAWING
				2:1	SHEET 3 OF 3