

BILL OF MATERIAL:

SUFFIX	LUER BODY	CONNECTOR	STEM
024	CLEAR POLYCARBONATE MAKROLON Rx1805 451118	CLEARPOLYCARBONATE MAKROLON Rx1805 451118	BLUE SILICON 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238
024R	CLEAR POLYCARBONATE MAKROLON Rx1805 451118 WITH RED COLORANT, NC31664830	CLEARPOLYCARBONATE MAKROLON Rx1805 451118	BLUE SILICON 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238
024FT	CLEAR POLYCARBONATE MAKROLON Rx1805 451118	CLEARPOLYCARBONATE MAKROLON Rx1805 451118	BLUE SILICON 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238
024RFT	CLEAR POLYCARBONATE MAKROLON Rx1805 451118 WITH RED COLORANT, NC31664830	CLEARPOLYCARBONATE MAKROLON Rx1805 451118	BLUE SILICON 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238