



DRAWING NO. C2455280XXX		REV B	Customer Drawing		REVISION HISTORY														
REV		DESCRIPTION (ECO)		DATE		APPROVED													
A		NEW DRAWING ECO 7199		9/23/16		KJS													
B		ADD ISO80369-7 TO NOTE, ECO 8017		4/28/18		AK													
<div>BILL OF MATERIAL:</div> <table><tr><td>SUFFIX</td><td>LUER BODY</td><td>CONNECTOR</td><td>STEM</td></tr><tr><td>24</td><td>CLEAR POLYCARBONATE MAKROLON Rx 1805 451118</td><td>CLEAR POLYCARBONATE MAKROLON Rx 1805 451118</td><td>BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238</td></tr><tr><td>24R</td><td>CLEAR POLYCARBONATE MAKROLON Rx 1805 451118 AND RED COLORANT</td><td>CLEAR POLYCARBONATE MAKROLON Rx 1805 451118</td><td>BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238</td></tr></table>								SUFFIX	LUER BODY	CONNECTOR	STEM	24	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238	24R	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118 AND RED COLORANT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238
SUFFIX	LUER BODY	CONNECTOR	STEM																
24	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238																
24R	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118 AND RED COLORANT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238																
MOLDED PARTS: GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA		UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		THIRD ANGLE PROJECTION		TITLE Tube End Port 2.8mm													
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION. ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.		DRAWN BY DMC		DATE 09/23/2016		HR® Halkey-Roberts®													
		APPROVED BY KJS		DATE		SIZE B													
		MATERIAL				DRAWING NO. C2455280XXX													
		FINISH				REV B													
		SCALE 4:1		DO NOT SCALE DRAWING		SHEET 2 OF 2													