


DRAWING NO. C2455205XXX		REV C	Customer Drawing		REVISION HISTORY			
REV	DESCRIPTION (ECO)			DATE	APPROVED			
B	ADD SUFFIX 24R ECO 7784			12/07/17	AK			
C	ADD ISO80369-7 TO NOTES ECO 8119			06/14/18	AK			

MATERIAL:

SUFFIX "XXX"	LUER BODY	CONNECTOR	STEM
24	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40
24R	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118, W/ RED COLORANT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40

DRAWN BY		DMC	DATE	09/01/2016
APPROVED BY			DATE	
MATERIAL				
FINISH				

MOLDED PARTS: GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA		UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005		THIRD ANGLE PROJECTION 	
PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION, ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.		SIZE B		TITLE HR® Halkey-Roberts® Tube End Port 2.0mm	
		SCALE 2:1		DRAWING NO. C2455205XXX SHEET 2 OF 2	