

DRAWING NO. C2455285XXX		REV C	Customer Drawing		REVISION HISTORY			
REV		DESCRIPTION (ECO)		DATE		APPROVED		
B		ADD 24R SUFFIX TO NOTE 1. ECO 7784		12/08/17		KJS		
C		ADD ISO80369-7 TO NOTES, ECO 8017		4/30/18		AK		

MATERIAL:

SUFFIX "XXX"	LUER BODY	CONNECTOR	STEM
24	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40
24R	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118 W/ RED COLORANT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40

DRAWN BY

DMC

DATE

9/26/16

APPROVED BY

DATE

MATERIAL

FINISH

MOLDED PARTS:

GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION. ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:

DIMENSIONS ARE IN INCHES


TOLERANCES:

ANGULAR: ± 1°

TWO PLACE DECIMAL ± .01

THREE PLACE DECIMAL ± .005

THIRD ANGLE PROJECTION



TITLE

HR® Halkey-Roberts®

SIZE

B

SCALE

2:1

DRAWING NO.

C2455285XXX

REV

C

SHEET

2 OF 2

DO NOT SCALE DRAWING