


BILL OF MATERIAL:

NO SUFFIX	CELCON M270
HF	HOSTAFORM MT24U01
GAM	(NATURAL POLYESTER) VALOX 325M
P	WHITE POLYPROPYLENE H0500HN
PBG,PG	GAMMA STABLE POLYPROPYLENE P5M6K-080

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY	M.S.	DATE		
APPROVED BY	J.W.	DATE	10/24/94	
MATERIAL				
FINISH				
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: $\pm 1^{\circ}$ TWO PLACE DECIMAL $\pm .01$ THREE PLACE DECIMAL $\pm .005$		THIRD ANGLE PROJECTION 		SIZE A
SCALE 2:1		HR[®] Halkey-Roberts[®] TUBING CLAMP		
		DRAWING NO. C340TCX		REV W
		DO NOT SCALE DRAWING		SHEET 2 OF 2