

DRAWING NO.
C340TCS2ABS



REV
A



Customer Drawing

Revision History

| REV | DESCRIPTION (ECO) | DATE | APPROVED |
|-----|---|---------|----------|
| A | UPDATE NOTES, REMOVE DISCONTINUED SUFFIX OPTIONS. NEW DRW. ECN 4057 | 5/21/09 | P.M. |
| | | | |



NOTES:


1. MATERIAL: ABS
2. CRITICAL DIMENSIONS ARE CIRCLED AND SUBJECT TO INSPECTION.
3. DIMENSIONS IN PARENTHESES ARE FOR REFERENCE ONLY.
4. ALL DIMENSIONS ARE IN INCHES WITH NOMINAL DIMENSIONS IN MILLIMETERS BELOW.
5. MAX TUBING DIA. (0.25")

QUALITY CRITERIA

O.D. FLASH 0.015 MAX. (0.4MM)

GATE VESTIGE 0.01 MAX. (0.3MM)

EMBEDDED PARTICLES, BUBBLES, SPOTS, STAINS: 0.3mm² MAX. NO MORE THAN 5 PER PART.

| | | | |
|------------------------------------|----------------|---|---------|
| DRAWN BY | DMC | DATE | 5/21/09 |
| APPROVED BY | - | DATE | - |
| MATERIAL | ABS#348, WHITE | | |
| FINISH | | | |
| UNLESS OTHERWISE SPECIFIED: | | THIRD ANGLE PROJECTION | |
| DIMENSIONS ARE IN INCHES | |  | |
| TOLERANCES: | | | |
| ANGULAR: ± 1° | | | |
| TWO PLACE DECIMAL ± .01 | | | |
| THREE PLACE DECIMAL ± .005 | | | |



Halkey-Roberts

TITLE **MEDI - CLAMP, SMALL**

| | | | | | |
|-------|----------|----------------------|--------------------|-------|----------|
| SIZE | A | DRAWING NO. | C340TCS2ABS | REV | A |
| SCALE | 2:1 | DO NOT SCALE DRAWING | | SHEET | 1 OF 1 |

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.