


BILL OF MATERIAL:

LUER BODY	Y CONNECTOR	STEM
CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	BLUE SILICONE 40 DUROMETER WACKER LR3003/40 DT COLOR, K75238

MOLDED PARTS:  
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL  
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY  
OF HALKEY-ROBERTS CORPORATION, ST. PETERSBURG, FL. ANY  
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN  
PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:  
DIMENSIONS ARE IN INCHES  
TOLERANCES:  
ANGULAR: ± 1°  
TWO PLACE DECIMAL ± .01  
THREE PLACE DECIMAL ± .005

THIRD ANGLE PROJECTION  


DRAWN BY  
DMC

APPROVED BY  
PM

MATERIAL

FINISH

DATE  
10/30/14

DATE

HR® Halkey-Roberts®

TITLE  
Y CONNECTOR SWABABLE

SIZE  
B

SCALE  
2:1

DRAWING NO.  
C2456240XX

DO NOT SCALE DRAWING

REV  
D

SHEET  
2 OF 2