



BILL OF MATERIAL:

1ST SUFFIX "X"	2ND SUFFIX "X"	LUER BODY	CONNECTOR	STEM
24	FT	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR POLYCARBONATE MAKROLON Rx 1805 451118	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40
50	FT	CLEAR EASTAR COPOLYESTER MN004	CLEAR EASTAR COPOLYESTER MN004	CLEAR SILICONE 40 DUROMETER WACKER LR3003/40