





NOTES:

1. LUER SYRINGE ACTIVATED TWO WAY CHECK VALVE.
COMPATIBLE WITH ISO 80369-7:2021.
2. HOUSING: WHITE POLYCARBONATE
STEM: SILICONE
2.5ml: LASER PRINT
SPRING: 302 STAINLESS STEEL PASSIVATE
3. VENDOR MUST SUPPLY CERTIFICATION THAT IDENTIFIES THE
CUSTOMER PART NUMBER, REVISION LEVEL, VENDOR LOT NUMBER,
IDENTIFICATION OF MATERIAL AND QUANTITY.
4. ALL COMPONENTS HAVE BEEN 100% LEAK TESTED TO A MAX
ALLOWABLE LEAKAGE RATE OF <0.055CC/MIN.
5. CRITICAL DIMENSIONS ARE CIRCLED AND SUBJECT TO INSPECTION.
6. DIMENSIONS IN PARENTHESES ARE FOR REFERENCE ONLY.
7. ALL DIMENSIONS ARE IN INCHES WITH NOMINAL DIMENSIONS IN
METRIC BELOW.



MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY
OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY
REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN
PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY	L.CHARBONNEAU	DATE	07/01/2025
APPROVED BY	J. SCHINIS	DATE	08/04/2025
MATERIAL	SEE NOTES		
FINISH			
UNLESS OTHERWISE SPECIFIED:		THIRD ANGLE PROJECTION	
DIMENSIONS ARE IN INCHES			
TOLERANCES:			
ANGULAR: ± 1°			
TWO PLACE DECIMAL ± .01			
THREE PLACE DECIMAL ± .005			

HR® Halkey-Roberts®			
TITLE			
CATHETER VALVE, 2.5ML			
SIZE	DRAWING NO.	REV	
A	CV2465420	A	
SCALE	DO NOT SCALE DRAWING		SHEET
5:1			1 OF 2

BILL OF MATERIAL:

LUER SLIP ADAPTER: WHITE POLYCARBONATE MAKROLON RX 1805013771
STEM: BLUE SILICONE
SPRING: 302 STAINLESS STEEL PASSIVATE
VALVE BODY: WHITE POLYCARBONATE MAKROLON RX 1805013771

MOLDED PARTS:
GATES, PARTING LINES, AND FLASH PER QUALITY CRITERIA

PROPRIETARY AND CONFIDENTIAL
THE INFORMATION CONTAINED IN THIS DRAWING IS SOLE PROPERTY OF HALKEY-ROBERTS CORPORATION ST. PETERSBURG, FL. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF HALKEY-ROBERTS CORPORATION IS PROHIBITED.

DRAWN BY L.CHARBONNEAU	DATE 07/01/2025	HR® Halkey-Roberts®	TITLE CATHETER VALVE, 2.5ML	SIZE A	DRAWING NO. CV2465420	REV A
APPROVED BY J. SCHINIS	DATE 08/04/2025					
MATERIAL SEE NOTES						
FINISH		THIRD ANGLE PROJECTION		SCALE 5:1		
UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ANGULAR: ± 1° TWO PLACE DECIMAL ± .01 THREE PLACE DECIMAL ± .005				DO NOT SCALE DRAWING		
				SHEET 2 OF 2		