



Nordson MEDICAL
801 West 71st Street
Loveland, CO 80538

+1.888.404.5837
+1.970.267.5200
nordsonmedical.com

October 29, 2019

Dear Valued Customer,

You are receiving this notice because your organization has purchased Value Plastics Inc. dba Nordson MEDICAL part numbers composed of Trinseo's Calibre Megarad 2081-15 FC03006 polycarbonate resin in the past 5 years. VP/Nordson MEDICAL part number and material suffix -9002.

We have been notified by the resin manufacture TRINSEO, that they have implemented changes to the resin. The mold release agent added to the resin has been changed by removing all animal derived components. There is no change in composition or CAS number associated with the animal free component change, and there should be no effect on the processing and performance of the finished products.

To assist our customers in this transition, Value Plastics Inc. dba Nordson MEDICAL has secured enough Calibre Megarad 2081-15 SC03006 polycarbonate resin to supply parts until the end of 2020, based upon previous and current demands. Implementation will be a running change with no scrap planned. Product manufactured with the new resin will be indicated on the Certificate of Conformance under the material information.

CURRENT CERTIFICATE OF CONFORMANCE IDENTIFIER OF RESIN GRADE:

Part #	Description	Customer Part #	Order Lot #	QTY
EXAMPLE Other Requirements:	Part Description	EXAMPLE	EXAMPLE	500
	Material Lot: EXAMPLE	Percent: 100	Material: TRINSEO Calibre Megarad 2081-15 FC030006	(9002)

FUTURE CERTIFICATE OF CONFORMANCE IDENTIFIER OF RESIN GRADE:

Part #	Description	Customer Part #	Order Lot #	QTY
EXAMPLE Other Requirements:	Part Description	EXAMPLE	EXAMPLE	500
	Material Lot: EXAMPLE	Percent: 100	Material: TRINSEO Calibre Megarad 2081-15 SC030006	(9002)

Supporting Documentation:

[EU MDR Statement](#)

[Technical Specification](#)

[MSDS](#)

[Regulatory Sheet](#)

[ISO 10993 Biocompatibility](#)

Thank you for your continued support and understanding as Nordson works to execute these changes implemented by our raw materials manufacturer. Value Plastics Inc. dba Nordson MEDICAL will assist our customers through this transition to the best of our ability and strives to ensure timely delivery of product.

Please contact your Account Manager, or Sales@nordsonmedical.com for additional assistance.

Sincerely,

Ben G. Johnson

Ben G. Johnson
North American Sales Manager