





Surface Area Data

Fluid Management



Part Series	Image	Part Numbers	Part Weight (Grams)	Wetted Surface Area (cm ²)
Straight and Straight Reducers		N220-6005	0.126	1.085
		N013-6005	0.216	1.732
		N035/013-6005	0.794	3.493
		N230-6005	0.227	1.988
		3060-6005	0.972	4.553
		N055-6005	1.724	7.091
		N055/013-6005	1.071	5.052
		N670-6005	2.495	12.969
		N070/013-6005	1.926	8.375
		N670/655-6005	1.950	9.542
		N680-6005	4.445	21.303
		N680/670-6005	3.778	17.155
		N680/655-6005	3.084	14.070
		N6110/6100-6005	15.858	55.936

Part Series	Image	Part Numbers	Part Weight (Grams)	Wetted Surface Area (cm ²)
Push in Plugs		PIP230-6005	0.136	0.663
		PIP60-6005	0.726	3.666
		PIP070-6005	1.629	4.940
Tee and Tee Reducers		T230-6005	0.311	2.795
		T60-6005	2.722	12.504
		T055-6005	1.568	9.498
		T670-6005	4.905	24.053
		T670/655-6005	4.653	17.875
		T680-6005	8.714	40.598
		T080/055-6005	8.800	27.968
		T680/655-6005	8.222	28.304
		T690-9	18.987	56.159
		T6100-6005	25.443	75.088
Cross Connectors		4PX655-6005	2.983	11.974
		4PX670-6005	8.208	28.399
		4PX680-6005	14.247	46.968
Y Connectors		Y230-6005	0.342	2.976
		Y630/655-6005	0.825	4.756
		Y260-6005	2.540	11.900
		Y670-6005	5.080	25.488
		Y680-6005	9.651	42.858
Luers		FTLL035-6005	0.454	2.639
		FTLL360-6005	0.606	3.953
Sanitary Fittings		SFMN670-6005	2.676	13.668
		SFMN680-6005	3.629	17.917
		SFMN6100-6005	6.336	22.513
		SFMX670-6005	7.218	23.493
		SFMX680-6005	7.872	28.474
		SFMX6100-6005	10.773	36.599
Elbow Connectors		SFMXC-6005	6.354	13.584
		L230-6005	0.224	1.867
		L60-6005	2.276	8.120
		L080-6005	7.738	22.980
		CF494-6005	18.558	67.514